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
The diagram shows a form titled "新生児スクリーニング採血ろ紙" (Newborn Screening Blood Filter Paper). The form is divided into several sections with callouts providing instructions:

- Header:** "1. 初回採血 2. 再採血 (回目)" (1. First blood draw 2. Re-blood draw (times)). Callout: "採血回数を記入してください 再検査の場合は2回目です" (Please enter the number of blood draws. In the case of a re-test, it is the 2nd time).
- Medical Institution:** "医療機関所在地" (Medical institution location) and "医療機関名" (Medical institution name). Callout: "市町村名まででも結構です" (It's fine to enter up to the city/town/village name).
- Residence:** "住所" (Address) with options for "1. 広島市()区" (1. Hiroshima City () Ward) and "2. 広島市以外の市町村()市・町・村" (2. City/Town/Village other than Hiroshima City () City/Town/Village). Callout: "病院内で区別が必要な場合は必ず科名、病棟/外来の区別をご記入ください" (If differentiation is needed within the hospital, please enter the department name, ward/outpatient clinic name).
- Parent Information:** "母フリガナ" (Mother's name in katakana) and "氏名" (Name). Callout: "カタカナでの報告になりますので、必ずフリガナをご記入ください 濁点などにもご注意ください" (The report will be in katakana, so please enter the name in katakana. Please be careful of punctuation marks, etc.). Example: "例)藤原→フジワラ/フジハラ 裕子→ユウコ/ヒロコ" (Example: Fujiwara → Fujiwara/Fujiwara Yoko → Yoko/Hiroko).
- Child Information:** "性別" (Sex) with "男・女" (Male/Female) and "未定" (Undetermined). Callout: "性別の記入をお願いいたします" (We request you to enter the sex). "在胎週数" (Gestational week) in "週" (Weeks). "出生日" (Date of birth) in "年 月 日" (Year, Month, Day). "出生時体重" (Birth weight) in "g" (grams). "哺乳開始日" (Date of breastfeeding start) in "年 月 日" (Year, Month, Day). "採血日" (Date of blood draw) in "年 月 日" (Year, Month, Day). "採血時体重" (Weight at blood draw) in "g" (grams). "抗生剤使用" (Antibiotic use) with "有・無" (Yes/No).
- Breastfeeding Status:** "哺乳" (Breastfeeding) with options "1. 良" (1. Good), "2. 不良" (2. Poor), "3. 殆ど哺乳不能" (3. Almost unable to breastfeed). Callout: "採血は生後5日目(日齢4日)以降、また哺乳開始から2日以上経過して採血してください" (Please draw blood after the 5th day of life (day 4) and after breastfeeding has started for more than 2 days).
- Results:** "結果" (Result) with "正常" (Normal) and "異常" (Abnormal). Callout: "記入例) 母(有)" (Example: Mother (Yes)).
- Footer:** "再採血必要" (Re-blood draw necessary), "精密検査必要" (Detailed examination necessary), "広島市医師会臨床検査センター" (Hiroshima City Medical Association Clinical Examination Center). Callout: "注目の場合はその旨を記載してください ただし、糖分のみ投与の場合は、『3. 殆ど哺乳不能』に丸をしてください" (Please note the reason if necessary. However, in the case of only sugar administration, please circle '3. Almost unable to breastfeed').

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